



Phantom Head Simulator Unit

Phantom head Simulator Unit includes

- Motorized Phantom Head with in-built Halogen light.
- Airotor control box with controls for Airotor.
- Micromotor (with 2 hand pieces) and 3-way syringe.
- Pneumatic Stool, Phantom head with half fibre glass body.
- Motorized up & down movement of body.
- Touch Pad control panel.
- Imported glass reflector.
- Typodont Jaw Set with 32 ivorine teeth.
- Enamel Finish.



For more Product range visit us at:

Bensons Surgico, Address : 42, Dlf Industrial Area, Kirti Nagar, New Delhi - 110 015 (INDIA)

Phone : +(91)-(11)-25932664/25158759, Tele Fax :+(91)-(11)-25932664/25464251

E-mail : bensons@bensonsmedical.com, info@bensons.in

Website : <http://www.bensonsmedical.com>